

PRISMA[®] MOBILE

PROGRESSIVE LENS

The ideal lens for users
of electronic devices.



Prisma Mobile

Improve your visual experience when using your **smartphone** or **tablet**.

Prisma Mobile is a progressive lens specifically designed for those who **frequently use electronic devices** such as cell phones, tablets, or smartwatches. Its layout provides wide visual fields for both near and distance vision, along with a soft, comfortable transition to the reading zone.



GREAT RESOLUTION THANKS TO DIGITAL RAY-PATH® TECHNOLOGY

Digital Ray-Path® is an innovative calculation technology that uses a design engine to compensate the lens with a simulation of the binocular eye-lens system. Every unique lens is individually calculated guaranteeing an adapted solution for any prescription and base curve.




SMART ADD TECHNOLOGY ALLOWS THE EYES TO BE RELAXED WHEN VIEWING SCREEN DISPLAYS

Prisma mobile II incorporates Smart Add, the technology that improves the patient's visual experience when using digital devices. This technology assists with changes in focus at different working distances in a more agile and accurate way when the patient is working or reading simultaneously with different screens.

SHORT CORRIDOR ALLOWS FOR QUICK TRANSITION TO NEAR VISION

Our constant use of smartphones and tablets has led to a greater need for switching between near and distance vision. Keeping the same corridor length, its unique progression profile offers a more comfortable transition from the distance zone to the reading zone.

PROGRESSION PROFILE



SMOOTH PROGRESSION PROFILE FOR GREATER COMFORT

The transition between distance and near vision is very smooth. Thanks to its unique layout, this lens offers a wide and comfortable reading zone.

TARGET & POSITIONING

- Premium personalized progressive lens for electronic device users
- Ideal for progressive lens wearers ages 40 and over, both experts and novices.



Digital lens



Personalized



Digital Ray-Path®



Smart Add



Multiple corridor



Short corridor available



Variable Inset



Wrap available

CONSUMER ADVANTAGES

With Prisma Mobile, wearers take advantage of these benefits:

REDUCED VISUAL FATIGUE

More relaxed eyes.

DYNAMIC VISION

Easy and agile switch from the near to intermediate zone.

IMPROVED FOCUS

Better ability to read text on digital displays.

NATURAL POSITION

Ergonomically comfortable.

FASTER ADAPTATION

The lens doesn't really need an adaptation period.

FULLY PERSONALIZED

Unique lenses for each patient.

HOW TO PRESCRIBE PRISMA MOBILE

PERSONALIZATION PARAMETERS

To provide 100% personalization, it is essential to include all personalization parameters unique to each wearer's prescription information. For orders that do not include personalization parameter data (i.e. pantoscopic angle, wrap angle, nasopupilar distance, back vertex distance and frame dimensions), the lens will be optimized by using default values.

MINIMUM FITTING HEIGHT:

14mm, 15mm, 16mm, 17mm, 18mm, 19mm, 20mm

MOUNTING INSTRUCTIONS

For proper positioning of the lens in the frame, the invisible engravings must be considered.

TRIANO Lab.
Building No 7766, st. 174 Second Industrial. City Riyadh. Saudi Arabia
Tel: +966 11 4986901- 011 4986905 - Fax: +966114777127
www.triano.com.sa